

| Scenario ID | JIRA ID name | Positive | Negative |
|-------------|---|--|--|
| 1 | Create a Personalized Planning Preference window that's displayable in Planning | <p>Pre-requisite: Ensure reference files are imported to CORI system.</p> <p>Log in as a surgeon-->Select SPP widget-->Select an implant, deformity and reference file--> Change preference values--> Note the values-->Quit--> Create a case--> Navigate to Combined planning (with & without gaps)--> Select SPP static widget--> Verify the following elements matches with SPP-Launcher screen:</p> <ol style="list-style-type: none"> 1. 3 static images of knee joint 2. selected deformity (non-editable) 3. Selected implant (non-editable) 4. Preference values set in SPP-Launcher (non-editable) 5. Reference drop down (Editable) 6. Reference values in orange color and displayed in closed brackets (non-editable) 7. Button to close the window 8. No Reset button | <p>Pre-requisite: Ensure reference files are imported to CORI system.</p> <p>Log in as a new surgeon-->Create a case--> Navigate to Combined planning (with & without gaps)--> Select SPP static widget--> Verify that the widget is disabled/ is hidden</p> |

| | |
|---|--|
| 2 | <p>Pre-requisite: Ensure reference files are imported to CORI system.</p> <p>Log in as a surgeon-->Select SPP widget-->Select an implant and deformity --> Change preference values--> Note the values-->Quit--> Create a case--> Navigate to Combined planning (with & without gaps)--> Select SPP static widget--> Verify the elements mentioned in cell G2 matches SPP-Launcher except the following:</p> <ol style="list-style-type: none"> 1. Reference files drop down displays "--choose--" and user can select any reference file 2. No reference values are displayed in closed brackets and in orange color. <p>Note: Detailed test cases for reference drop down is covered in RIOSI-929 so not covering in this sheet.</p> |
| 3 | <p>Execute G2 test scenario for surgeon created in 1.7 or earlier.</p> <p>Pre-requisite: Execute G2 TC</p> <p>Select SPP widget--> Change implant, deformity and reference files--> Quit--> Resume the same case--> Navigate to Combine dPlanning--> Verify the elements mentioned in G2 matches with SPP Launcher</p> |

| | | | |
|----|--|--|--|
| 4 | | <p>Pre-requisite: Ensure reference files are imported to CORI system.</p> <p>Log in as a surgeon-->Select SPP widget-->Select only implant--> Quit--> Create a case--> Navigate to Combined planning (with & without gaps)--> Select SPP static widget--> Verify the following elements matches with SPP-Launcher screen:</p> <ol style="list-style-type: none"> 1. 3 static images of knee joint 2. selected deformity displays "--Choose--" (non-editable) 3. Selected implant matches SPP Launcher (non-editable) 4. Preference values set in SPP-Launcher (non-editable and displays implant default) 5. Reference drop down (Editable) 6. No Reference values in orange color are displayed in closed brackets (non-editable) 7. Button to close the window 8. No Reset button <p>Note: Execute same TC for deformity</p> | |
| 5 | | | |
| 6 | | | |
| 7 | | | |
| 8 | | | |
| 9 | | | |
| 10 | | | |
| 11 | | | |
| 12 | | | |
| 13 | | | |
| 14 | | | |
| 15 | | | |
| 16 | | | |
| 17 | | | |
| 18 | | | |
| 19 | | | |
| 20 | | | |
| 21 | | | |
| 22 | | | |
| 23 | | | |
| 24 | | | |
| 25 | | | |
| 26 | | | |
| 27 | | | |
| 28 | | | |
| 29 | | | |

| | | | |
|----|--|--|--|
| 30 | | | |
| 31 | | | |
| 32 | | | |
| 33 | | | |
| 34 | | | |
| 35 | | | |
| 36 | | | |
| 37 | | | |
| 38 | | | |
| 39 | | | |
| 40 | | | |
| 41 | | | |
| 42 | | | |
| 43 | | | |
| 45 | | | |
| 46 | | | |
| 47 | | | |
| 48 | | | |
| 49 | | | |
| 50 | | | |
| 51 | | | |
| 52 | | | |

| Boundary Analysis | State Transition | Combination Testing |
|-------------------|---|---|
| NA | <p>Pre-requisite: Execute G2 TC</p> <p>Navigate to IRF screen via Admin Login--> Remove reference files-->Navigate to surgeon login-->Select SPP widget--> Change implant, deformity--> Quit--> Resume the same case--> Navigate to Combine dPlanning--> Verify the elements mentioned in G2 matches with SPP Launcher</p> | <p>Pre-requisite: Ensure reference files are imported to CORI system.</p> <p>Log in as a surgeon-->Select SPP widget-->Select an implant, deformity and reference file--> Change preference values--> Note the values-->Quit--> Create a case--> Navigate to Combined planning (with & without gaps)--> Select SPP static widget--> Verify the elements mentioned in G2 matches SPP Launcher screen--> Quit case--> Navigate to SPP Launcher screen --> Change all values to "--choose--"--> Quit--> Resume case--> Navigate to Combined Planning (with and without gaps)--> Verify that SPP widget is disabled/hidden</p> |

Pre-requisite: Ensure reference files are imported to CORI system.

Log in as a surgeon-->Select **SPP widget**-->Select an **implant and deformity** --> **Change** preference values--> Note the values-->Quit--> Go back to SPP screen--> Change all values to "--choose"--> Quit--> Create a TKA case--> Navigate to combined planning (with and without gaps)--> Select SPP static widget--> Veify that SPP widget is disabled

Pre-requisite: Ensure reference files are imported to CORI system.

Log in as a surgeon-->Select **SPP widget**-->Select an **implant and deformity** --> **Change** preference values--> Note the values-->Quit--> **Go back** to SPP screen--> **Change all** values --> **Quit**--> **Create** a TKA case--> Navigate to **combined planning (with and without gaps)**--> Select SPP static widget--> Verify that elements mentioned in G2 matches the latest changes in SPP Launcher

[illegible]

[illegible]

| Remarks | Linked JIRA ID |
|---------|--|
| | Linked Epic: https://smith-nephew.atlassian.net/browse/RIOSI-930 |

| |
|--|
| |
| |

| |
|--|
| |
|--|

[illegible]

[illegible]