

Scenario ID	JIRA ID name	SW Requirement	SY Requirement	Positive
1	Add an option to reset to preference default	NA(No requirement present)	NA(No requirement present)	<p>Pre-requisite: Ensure to change implant [Journey II XR], deformity[<7 Varus, <10 flex] and change preference values in SPP screen.</p> <p>Create a TKA case with Journey II XR and simulate deformity of [<7 varus, <10 flex]-->Navigate till Combined Planning (with/without gaps)--> Select more option widget--> Verify that there are two options for reset, Reset to personalized femur plan and Reset to personalized tibia plan-->Verify that they are enabled-->Verify that the button design matches other buttons present in More options widget</p>

2	NA(No requirement present)	NA(No requirement present)	Continue from G2 --> Navigate back to Combined planning screen-->Record resection depth, varus/valgus values, flexion/extension and slope values from all implant views--> Change implant position for only one parameter each for both femur and tibia --> Select more options --> Select to Reset personalized femur plan --> Verify that the resection depth, varus/valgus values return back to preferences set in SPP screen -->Navigate to Plan Tibia using widget --> Verify that the resection depth and slope values displayed are same as preference values set in SPP screen
---	----------------------------	----------------------------	--

3	NA(No requirement present)	NA(No requirement present)	<p>Pre-requisite: Ensure to change implant [Journey II XR], deformity[<7 Varus, <10 flex] and record preference parameters in SPP screen.</p> <p>Create a TKA case with Journey II XR--> Simulate deformity [<7 varus, <10 flex]--> Navigate to Combined planning screen--> Change the implant position--> Select more widgets--> Select Reset to personalized femur plan--> Verify that the resection depth values, varus/valgus and flexion extension values for femur are reset to the preference values saved in SPP--> Select more widgets--> Select Reset Femur Position button--> Verify that the resection depth values, varus/valgus and flexion extension values for femur are reset to implant default which is different from preference values</p>
4	NA(No requirement present)	NA(No requirement present)	<p>Continue from G4-->Select more widgets-->Select Reset to personalized tibia plan--> Verify that the resection depth values, varus/valgus and slope values for tibia are reset to the preference values saved in SPP--> Select more widgets--> Select Reset Tibia Position button--> Verify that the resection depth values, varus/valgus and slope values for tibia are reset to implant default which is different from preference values</p>

5	NA(No requirement present)	NA(No requirement present)	Continue from G5 --> Navigate to Plan Tibia through more options widget--> Change implant positions--> Select Reset Tibia position --> Navigate to Combined Planning screen --> Select more widgets -->Select Reset to personalized tibia plan --> Verify that the resection depth values, varus/valgus and slope values for tibia are reset to the preference values saved in SPP
6	NA(No requirement present)	NA(No requirement present)	Continue from G6 -->Navigate to Plan Femur through more options widget --> Change implant positions-->Navigate to Combined Planning screen -->Select more options widget --> Select Reset Femur position -->Verify that the implant positions change to implant defaults and does not match with preference values saved in SPP-->Select more options widget --> Select Reset to personalized femur plan --> Verify that resection depth values, varus/valgus values, flexion extention values changes to preference values set in SPP

7	NA(No requirement present)	NA(No requirement present)	<p>Pre-requisite: Ensure to create a TKA case with Legion PS and deformity [> 7 Valgus > 10 flex] and complete it in 1.7 software version or earlier</p> <p>Uparde software version to CORI 2.0--> Navigate to SPP--> Select Legion PS and deformity [> 7 valgus < 10 flex]--> Change parameters--> Record the preference parameters--> Quit from SPP--> Resume case--> Navigate to Combined Planning screen--> Select more options widget--> Verify that Reset to personalized femur and tibia buttons are disabled--> Navigate back and clear all points till Robotic drill connection state--> Simulate deformity [> 7 valgus < 10 flex] --> Navigate till combined planning screen--> Select more options widget--> Verify that Reset to personalized femur and tibia buttons are enabled--> Change implant positions for both femur and tibia--> Select more options widget--> Click Reset to personalized femur and tibia buttons--> Verify that the implant positions values matches with preference values saved in SPP</p>
8			
9			
10			
11			
12			
13			

14				
15				
16				
17				
18				
19				
20				
21				
22				
23				
24				
25				
26				
27				
28				
28				
29				
30				
31				
32				
33				
34				
35				
36				
37				
38				
39				
40				
41				
42				
43				
45				
46				
47				
48				
49				
50				
51				
52				

Negative	Boundary Analysis
<p>Pre-requisite: Navigate to SPP and select implant [Journey II CR]--> Quit</p> <p>Create a TKA case with Journey II CR-->Navigate to Combined Planning--> Select more options widget--> Verify that Reset to personalized femur plan and tibia plan buttons are disabled--> Change implants other than Journey II CR and check in more widgets that Reset to personalized femur plan and tibia plan buttons are disabled</p>	NA

Pre-requisite: Navigate to SPP and select deformity [<7 varus, < 10 flex]--> Quit

Create a TKA case with Anthem PS-->Simulate deformity [<7 varus, < 10 flex]-->Navigate to **Combined Planning**--> Select **more options** widget--> Verify that Reset to personalized femur plan and tibia plan buttons **are disabled**--> **Change implants** other than Anthem PS and check in **more widgets that Reset to personalized femur plan and tibia plan buttons are disabled**

<p>Pre-requisite: Ensure to create a new surgeon and don't access SPP Launcher</p> <p>Log in as surgeon--> Create a new TKA case with Anthem PS--> Navigate to Combined Planning--> Select more options widget--> Verify that Reset to personalized femur plan and tibia plan buttons are disabled--> Change implants in combined planning screen --> Verify for each implant that Reset to personalized femur plan and tibia plan buttons are disabled</p>	
<p>Continue from H3--> quit from case--> Navigate to SPP screen--> Change implant to Legion PS, deformity [<7 varus, >10 flex]--> Change preference parameters--> Quit--> Reset case--> Simulate deformity [<7 varus, >10 flex]--> Navigate to Combined planning screen--> Select more options widget--> Verify for each implant that Reset to personalized femur plan and tibia plan buttons are disabled--> Change implants in Combined planning screen and verify that except for Legion PS the reset buttons for personlaized femur and tibia plan are disabled</p>	

<p>Continue from G3--> Navigate to Plan Femur--> Verify that there is only one Reset button present for resetting to implant default--> Navigate to Plan Tibia through Combined Planning-->Verify that there is only one Reset button present for resetting to implant default</p>	
<p>Pre-requisite: Ensure to change implant [Journey II BCS], deformity[<7 Varus, >10 flex] and change preference values in SPP screen.</p> <p>Log in as surgeon--> Create a new TKA case with Journey II BCS--> Simulate deformity of [> 7 valgus, >10 flex]-->Navigate to Combined Planning--> Select more options widget--> Verify that Reset to personalized femur plan and tibia plan buttons are disabled</p>	


<p>Pre-requisite: Ensure to change implant [Journey II CR], deformity [<7 Varus, >10 flex], Change preference paraameters. Change implant [Legion CR], deformity [>7 Valgus, > 10 flex], change preference parameters.</p> <p>Start TKA case with Legion CR--> Simulate deformity [<7 Varus, >10 flex]--> Navigate to combined planning-->Select more options widget--> Verify that Reset to personalized femur plan and tibia plan buttons are disabled--> Change implant on Combined planning screen to Journey II CR-->Select more options widget--> Verify that Reset to personalized femur plan and tibia plan buttons are enabled--> Navigate back till Robotic drill connection state and clear points--> Simulate deformity [>7 Valgus, > 10 flex]--> Navigate to combined planning-->Select more options widget--> Verify that Reset to personalized femur plan and tibia plan buttons are disabled</p>	
<p>Pre-requisite: Ensure to create a new surgeon and start the following test steps without accessing SPP widget.</p> <p>Start TKA case with Legion CR--> Simulate deformity [<7 Varus, >10 flex]--> Navigate to combined planning-->Verify that Reset to personalized femur plan and tibia plan buttons are disabled--> Click on the buttons and verify that no change is observed on implant positions for both femur and tibia.</p>	

[illegible]

State Transition	Combination Testing
<p>Continue from H5-->Navigate back and clear all points till Robotic Drill connection state--> Collect good points and simulate deformity [<7 Varus, >10 flex] --> Navigate to Combined Planning screen--> Select more options widget--> Verify that Reset to personalized femur plan and tibia plan buttons are enabled</p>	<p>Pre-requisite: Ensure to change implant [Anthem PS], deformity[<7 Varus, <10 flex] and change preference values in SPP screen.</p> <p>Create a TKA case with Anthem PS and simulate deformity of [<7 varus, <10 flex]-->Navigate till Combined Planning (with/without gaps)--> Chnage implant positions-->Select more option widget--> Verify that there are two options for reset, Reset to personalized femur plan and Reset to personalized tibia plan-->Verify that they are enabled-->Verify that the button design matches other buttons present in More options widget-->Change implant positions-->Click on Reset to personalized plan buttons for femur and tibia-->Verify that the implant positions reset to preference values set in SPP -->Click Quit from case--> Navigate to SPP--> change preference paramters--> Quit from SPP--> Resume case--> Navigate to Combined planning screen--> Select more option widget-->Click on Reset to personalized plan buttons for femur and tibia-->Verify that the implant positions reset to preference values set in SPP</p>

Pre-requisite: Ensure to change implant [Genesis II PS], deformity[>7 Varus, >10 flex] and change preference values in SPP screen.


Create a **TKA case with Genesis II PS and simulate deformity of [>7 varus, >10 flex]**-->Navigate till **Combined Planning (with/without gaps)**--> Change implant positions--> Change implant positions-->Select more option widget-->Click on **Reset to personalized plan buttons for femur and tibia**-->Verify that the implant positions reset to preference values set in SPP -->**Click Quit** from case--> Navigate to SPP--> Reset the preference plan--> Record preference values-->**Quit** from SPP--> **Resume** case--> Navigate to Combined planning screen--> Select more option widget-->Click on Reset to personalized plan buttons for femur and tibia-->Verify that the implant positions reset to preference values set in SPP



[illegible]

[illegible]

Remarks	Linked JIRA ID
	Linked Epic: https://smith-nephew.atlassian.net/browse/RIOSI-933



--

--

[illegible]

