

Scenario ID	JIRA ID name	Scenario Description	Feature
1	Update back button boundary shape in all states	Test Description: TBD Acceptance Criteria: TBD Testing Strategy: TBD Tools Used: N/A NOTE: 1. All the 3 monitors such as CORI monitor, Tablet and display monitors are connected to the console during execution. 2. Refer Appendix sheet to verify the scenarios. 3. Ensure the language is chosen by default i.e., English. 4. Ensure the PHI is set to default i.e., PHI access enabled. 5. Refer State Transition sheet for the state transition states and the list of states to be verified.	UI/UX Update
2			
3			
4			
5			
6			
7			
8			
9			

10			
----	--	--	--

Scenario Type	Scenario
Positive	Navigate to all CORI states via Surgeon Login --> Verify "Back" button is updated (shape, size and position), the button is functional (by touch - single press) and feasible to touch (long press and swipe) with 3 monitors.
Positive	Navigate to all CORI states via Surgeon Login --> Verify "Back" button is updated (shape, size and position), the button is functional (by foot pedal press - single press) and feasible to press by foot pedal (long press and multiple press) with 3 monitors.
Negative	Navigate to all CORI states via Surgeon Login --> Verify "Back" button is updated (shape, size and position), the button is non-functional (by touch - single press) and not feasible to touch (long press and swipe) with 3 monitors when notification occurs and the "Back" button is disabled.
Negative	Navigate to all CORI states via Surgeon Login --> Verify "Back" button is updated (shape, size and position), the button is non-functional (by foot pedal press - single press) and not feasible to press by foot pedal (long press and multiple press) with 3 monitors when notification occurs and the "Back" button is disabled.
State Transition	Navigate to all CORI states via Surgeon Login by foot pedal press --> Verify "Back" button is updated (shape, size and position) --> Navigate to all the states from the last state (say post-op gap assessment state) by touchscreen only (reverse workflow) --> Verify "Back" button is updated (shape, size and position) & the button is functional (by touch - single press) with 3 monitors.
State Transition	Navigate to all CORI states via Surgeon Login by touchscreen --> Verify "Back" button is updated (shape, size and position) --> Navigate to all the states from the last state (say post-op gap assessment state) by foot pedal press only (reverse workflow) --> Verify "Back" button is updated (shape, size and position) & the button is functional (by foot pedal press - single press) with 3 monitors.
State Transition	Navigate to all CORI states via Surgeon Login by foot pedal press --> Verify "Back" button is updated (shape, size and position) --> From the last state (say post-op gap assessment state) navigate to each state from all the states mentioned in state transition table via widget & Robotic Drill Disconnection --> Verify "Back" button is updated (shape, size and position) & the button is functional (by foot pedal press - single press) with 3 monitors.
State Transition	Navigate to all CORI states via Surgeon Login by touchscreen --> Verify "Back" button is updated (shape, size and position) --> From the last state (say post-op gap assessment state) navigate to each state from all the states mentioned in state transition table via widget & Robotic Drill Disconnection --> Verify "Back" button is updated (shape, size and position) & the button is functional (by touch - single press) with 3 monitors.
Combination Testing	Resume the case created at Scenario 1 --> Navigate to all CORI states after resume operation via foot pedal press --> Verify "Back" button is updated (shape, size and position) & the button is functional (by foot pedal press - single press) with 3 monitors --> Navigate to all states using back button via touchscreen (reverse workflow) --> Verify "Back" button is updated (shape, size and position) & the button is functional (by touch - single press) with 3 monitors.

Combination Testing	Restart the case created at Scenario 1 --> Navigate to all CORI states after restart the operation via foot pedal press --> Verify "Back" button is updated (shape, size and position) & the button is functional (by foot pedal press - single press) with 3 monitors --> Navigate to all states using back button via touchscreen (reverse workflow) --> Verify "Back" button is updated (shape, size and position) & the button is functional (by touch - single press) with 3 monitors.
------------------------	---

PGH Feedback	Linked JIRA ID
	Linked Epic: https://smith-nephew.atlassian.net/browse/CORI-502



