

Scenario ID	JIRA ID name	Scenario Description	Feature
1	Update exit icon in all states	Test Description: TBD Acceptance Criteria: TBD Testing Strategy: TBD Tools Used: N/A NOTE: 1. All the 3 monitors such as CORI monitor, Tablet and display monitors are connected to the console during execution. 2. Refer Appendix sheet to verify the scenarios. 3. After exiting or quitting the case, resume the same case to check for other states in the workflow for the same scenario. 4. "Exit" or "Quit" icon is always enabled. It is masked when notification pops up. 5. Ensure the language is chosen by default i.e., English. 6. Ensure the PHI is set to default i.e., PHI access enabled. 7. Refer State Transition sheet for the state transition states and the list of states to be verified.	UI/UX Update
2			
3			
4			
5			
6			
7			
8			

Scenario Type	Scenario
Positive	Navigate to all CORI states via Surgeon Login --> Verify "Exit" icon is updated (shape, size and position), the button is functional (by touch - single press) & feasible to touch (long press and swipe) with 3 monitors.
Negative	Navigate to all CORI states via Surgeon Login --> Verify "Exit" icon is visible & updated (shape, size and position), the button is non-functional (by touch - single press) & not feasible to touch (long press and swipe) with 3 monitors when notification occurs and "Exit" icon is disabled.
State Transition	Navigate to all CORI states via Surgeon Login by foot pedal press --> Verify "Exit" icon is updated (shape, size and position) & the button is functional (by touch - single press) with 3 monitors --> From the last state (say post-op gap assessment state) navigate using Back button to all states by touchscreen only (reverse workflow) --> Verify "Exit" icon is updated (shape, size and position) & the button is functional (by touch - single press) with 3 monitors.
State Transition	Navigate to all CORI states via Surgeon Login by touchscreen --> Verify "Exit" icon is updated (shape, size and position) & the button is functional (by touch - single press) with 3 monitors --> From the last state (say post-op gap assessment state) navigate using Back button to all states by foot pedal press only (reverse workflow) --> Verify "Exit" icon is updated (shape, size and position) & the button is functional (by touch - single press) with 3 monitors.
State Transition	Navigate to all CORI states via Surgeon Login by foot pedal press --> Verify "Exit" icon is updated (shape, size and position) --> From the last state (say post-op gap assessment state) navigate to each state from all the states mentioned in state transition table via widget & Robotic Drill Disconnection --> Verify "Exit" icon is updated (shape, size and position) & the button is functional (by touch - single press) with 3 monitors.
State Transition	Navigate to all CORI states via Surgeon Login by touchscreen --> Verify "Exit" icon is updated (shape, size and position) --> From the last state (say post-op gap assessment state) navigate to each state from all the states mentioned in state transition table via widget & Robotic Drill Disconnection --> Verify "Exit" icon is updated (shape, size and position) & the button is functional (by touch - single press) with 3 monitors.
Combination Testing	Resume the case created at Scenario 1 --> Navigate to all CORI states after resume operation via foot pedal press --> Verify "Exit" icon is updated (shape, size and position) & the button is functional (by touch - single press) with 3 monitors --> Navigate to all states using back button via touchscreen (reverse workflow) -->Verify "Exit" icon is updated (shape, size and position) with 3 monitors.
Combination Testing	Restart the case created at Scenario 1 --> Navigate to all CORI states after restart the operation via foot pedal press --> Verify "Exit" icon is updated (shape, size and position) & the button is functional (by touch - single press) with 3 monitors --> Navigate to all states using back button via touchscreen (reverse workflow) -->Verify "Exit" icon is updated (shape, size and position) with 3 monitors.

PGH Feedback

[illegible]

Linked JIRA ID
<p>Linked Epic:</p> <p>https://smith-nephew.atlassian.net/browse/CORI-502</p>