

Scenario ID	JIRA ID name	Scenario Description	Feature
1	Redesign icons	<p> Test Description: TBD Acceptance Criteria: TBD Testing Strategy: TBD Tools Used: N/A NOTE: 1. All the 3 monitors such as CORI monitor, Tablet and display monitors are connected to the console during execution. 2. Refer Appendix sheet to verify the scenarios. 3. Ensure the language is chosen by default i.e., English. 4. Ensure the PHI is set to default i.e., PHI access enabled. 5. Refer State Transition sheet for the state transition states, the list of redesigned icons and the list of states to be verified. 6. Refer Redesign Icons sheet to see the redesigned icon image. </p>	UI/UX Update
2			
3			
4			
5			
6			
7			
8			

Scenario Type	Scenario
Positive	Navigate to all CORI states via Surgeon Login --> Verify all the icons are redesigned (shape, size and position), the button is functional (by touch - single press) & feasible to touch (long press and swipe) with 3 monitors.
Negative	Navigate to all CORI states via Surgeon Login --> Verify all the icons are visible & redesigned (shape, size and position) & the button is non-functional (by touch - single press) with 3 monitors when notification occurs and all the icons are disabled.
State Transition	Navigate to all CORI states via Surgeon Login by foot pedal press --> Verify all the icons are redesigned (shape, size and position) & the button is functional (by touch - single press) with 3 monitors --> From the last state (say post-op gap assessment state) navigate using Back button to all states by touchscreen only (reverse workflow) --> Verify all the icons are redesigned (shape, size and position) & the button is functional (by touch - single press) with 3 monitors.
State Transition	Navigate to all CORI states via Surgeon Login by touchscreen --> Verify all the icons are redesigned (shape, size and position) & the button is functional (by touch - single press) with 3 monitors --> From the last state (say post-op gap assessment state) navigate using Back button to all states by foot pedal press only (reverse workflow) --> Verify all the icons are redesigned (shape, size and position) & the button is functional (by touch - single press) with 3 monitors.
State Transition	Navigate to all CORI states via Surgeon Login by foot pedal press --> Verify all the icons are redesigned (shape, size and position) --> From the last state (say post-op gap assessment state) navigate to each state from all the states mentioned in state transition table via widget & Robotic Drill disconnection --> Verify all the icons are redesigned (shape, size and position) & the button is functional (by touch - single press) with 3 monitors.
State Transition	Navigate to all CORI states via Surgeon Login by touchscreen --> Verify all the icons are redesigned (shape, size and position) --> From the last state (say post-op gap assessment state) navigate to each state from all the states mentioned in state transition table via widget & Robotic Drill disconnection --> Verify all the icons are redesigned (shape, size and position) & the button is functional (by touch - single press) with 3 monitors.
Combination Testing	Resume the case created at Scenario 1 --> Navigate to all CORI states after resume operation via foot pedal press --> Verify all the icons are redesigned (shape, size and position) & the button is functional (by touch - single press) with 3 monitors --> Navigate to all CORI states using back button via touchscreen (reverse workflow) --> Verify all the icons are redesigned (shape, size and position) with 3 monitors.
Combination Testing	Restart the case created at Scenario 1 --> Navigate to all CORI states after restart the operation via foot pedal press --> Verify all the icons are redesigned (shape, size and position) & the button is functional (by touch - single press) with 3 monitors --> Navigate to all CORI states using back button via touchscreen (reverse workflow) --> Verify all the icons are redesigned (shape, size and position) with 3 monitors.

PGH Feedback

Linked JIRA ID
<p>Linked Epic:</p> <p>https://smith-nephew.atlassian.net/browse/CORI-502</p>